

3Lit Mini Sterilizer

Compact, Speedy, Big capacity

Why is X3 compact?



The configuration is scientific and ratioal

- Internal reservior / draining tank
- Build-in printer / SD menory




Why is X3 speedy?



Top performance of complete Class B sterilization

Cycle with drying is 13 minutes.(datas from our lab)



Why is capacity big?



The chamber is designed specially of rectangle shape.



1 implant case




2 trays

2 trays or 1 implant case can be accomodated



More than those:





Side-open door, good sealability and high safety



X3 mini Technical Data					
Power supply	230V±23V, 50Hz,2900VA				
Rack & trays	Standard 2 trays				
External dim.WxHxD	287 x 248 x595 mm				
Chamber dim. WxHxD	175 x 55 x 280 mm				
N.Weight/G.Weight	24kgs/29kgs				
Water tank	1.0 Lit				
Waste tank	0.8Lit				

Cycles	Sterilization temperature (°C)	Holiding time(min)	Pre-vacuum times	Drying Time(min)	Remarks	
134°C Jniversal Quick	134	4	3	4	Preset	
34°C Universal	134	4	3	6	Preset	
34°C Solid	134	4	1	2		
21°C Universal	121	4	3	8	Preset	
3-D/Helix test	134	3.5	3	4	Test program	
/acuum test		Vacuum to -80kPa,holding 15min			Test program	
Drying				5	Preset drying	
						*